HUNTLEIGH

dopplex[®] DR4

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Mode d'emploi

Bruksanvisning

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anwijzing

INSTRUCTIONS FOR USE

Bruksanv

matları

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Version 4

ENGLISH

Contents	Page No.
1. Introduction	3
1.1 Product Overview	3
1.2 Package Contents	3
1.3 Hardware Requirements	4
1.4 Installation of Software	4
2. Getting Started	6
2.1 Starting the Program	6
2.2 Patient Entry Information	6
2.3 Configure Measurement	7
2.3.1 Brachial Pressure	8
2.4 Doppler Communications	9
2.5 Doppler Application	10
2.6 Doppler Calculations	11
2.6.1 Auto Gain Control	12
2.7 VPPG Application	12
2.8 ABI Application	14

SERVICE HELPLINE

In Europe: +44 (0) 2920 485885
 email: cardiff.service@huntleigh-diagnostics.co.uk

In USA: 1800.323.1245 (option 2)

dopplex, Huntleigh and 'H' logo are registered trademarks of Huntleigh Technology Ltd.

Microsoft and Windows are trademarks of Microsoft Corporation

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

1. Introduction

1.1 Product Overview

Thank-you for purchasing the Huntleigh Healthcare dopplex® DR4 software package.

The DR4 has been designed to operate with the dopplex® MD2, RD2 and MD200 bi-directional Dopplers, and the dopplex® Ability. This package will run on a Pentium computer and will display bi-directional velocity/time waveforms obtained from various vessels. It will also display pressures, ABIs and PVR waveforms generated from the dopplex Ability. Data can be stored on a disk and a printout can be obtained for patient notes.

The VPPG application can record and store VPPG waveforms when connected to the RD2.

Doppler parameters can be calculated on new and existing blood flow waveforms.

The software requires Microsoft™ Windows™ XP, Vista, W7 or W8 to be installed on the computer prior to use. (See 'Installation of Software' section).

An extensive Help facility is available to explain these additions and how to operate the software.

1.2 Package Contents

The following items should be contained within the software package:

1. License agreement
2. CD ROM containing the software, Instructions for Use and USB Adaptor Installation Guide.
3. RS232 serial interface cable - for connection to the MD2/RD2/MD200
4. USB Serial Adaptor with driver CD

(Note: a USB cable is provided with dopplex Ability)

1.3 Hardware Requirements

The minimum hardware requirement is a Pentium PIII computer with 256MB of memory. The complete program will require at least 100MB of hard disk space.

1.4 Installation of Software

Microsoft Windows must be installed on your computer prior to installing DR4. If you require a copy of Windows, this can be obtained from your local computer dealer. To install Windows, follow the installation guide provided by Microsoft.

To install DR4 onto your hard disk:

1. Insert the CD ROM into your CD Rom drive. The program will automatically start Installation.
2. As the installation proceeds, DR4 prompts you to confirm or to provide information. The installation program provides default information, which you can accept or change.
3. The program will then be loaded onto the hard disk of your computer and the message '**Setup Completed**' will appear.
4. Any previous version will not be overwritten by DR4.
5. Existing patient files recorded with DR3 will be automatically moved to the 'DR4data' file and can then be viewed with DR4.
DR2 and DR3 patient files are fully compatible with DR4.
6. If you encounter any difficulty with installation of the software, telephone our Helpline number shown on the Contents page or contact your supplier.
7. The software has been designed to be used either with or without a mouse connected.

8. Connect the cable provided to the 9way RS232 serial port of your computer and to the **MD2/RD2/MD200**. If your computer does not have a 9 way serial port, use the USB-Serial adaptor). Follow the Installation instructions supplied on the CD.
9. To obtain maximum benefit from this software package, the entire Help feature should be read carefully.

2. Getting Started

2.1 Starting the Program

To start the DR4 program, double click on the  icon , or click once and press <Enter>.

The introductory window will be displayed as shown below:-

2.2 Patient Information Entry

Select the  icon and the '*Patient Information*' window will appear.

If the software is being run for the first time after installation, a new set-up file will be created.

If the MD2/RD2/MD200 lead is connected to an RS232 port that is not available for use by Windows, the following message will appear:

The 'Patient Information' dialog box contains the following fields and controls:

- Surname: [Text Box]
- Other Names: [Text Box]
- Street: [Text Box]
- Town: [Text Box]
- I.D.: [Text Box]
- Insurance Number: [Text Box]
- Height: [Text Box]
- Weight: [Text Box]
- Date of Birth: [Text Box]
- Gender: Male Female
- Test Date: 11/06/04
- GP: [Dropdown Menu]
- Reimbursement: [Dropdown Menu]
- Referred by: [Dropdown Menu]
- Performed by: [Dropdown Menu]
- Hospital: [Text Box]
- Location: [Text Box]
- Existing Conditions: [Text Area]
- Doctors Summary: [Text Area]
- Buttons: OK, Cancel, Help

Enter the information in the boxes, using the **<Tab>** key or mouse to move from box to box. Note that the **'Hospital'** and **'Location'** boxes can only be filled in from the **'Setup Menu - Location'**. Click on OK, or press **<Enter>** to accept the data.

2.3 Configure Measurement

Select the  icon and the Configure Measurement window will appear.

Select the type of display and the required settings, and click **OK** or press **<Enter>**.

The 'Configure Measurements' dialog box includes the following sections:

- Display:** Leg Arteries (dropdown)
- Preset:** API, Seg Pressure, Arterial (checkboxes)
- Waveforms:** Waveforms, Pressures, Both (radio buttons)
- Brachial:** Brachial (checkbox)
- Measurement Selection:** 12 numbered checkboxes (1-12) surrounding a leg diagram.
- Measurements Order:** a, b, c, d, e (radio buttons) and a list of order sequences.
- Cal. Togs:** Cal. Togs, Combined Waveforms, Invert (checkboxes)
- Auto Box Height:** Auto Box Height (checkbox)
- Default trace time (s):** 2.5, 5, 10, 20 (radio buttons)
- Default gain:** 1, 2, 4, 8, > (radio buttons)
- Buttons: OK, Cancel, Help


2.3.1 Brachial Pressure

To enter the Brachial Systolic Pressure, select the  icon.

Measure the patient's systolic blood pressure and enter the reading (in mmHg) in the window shown below.

Click on *OK* once again or press *<Enter>*.


Switch on the MD2/RD2/MD200 and ensure that *'DOPPLER ON'* is indicated at the bottom of the screen. If *'No Doppler'* is present, follow the instructions in the *'Doppler Communications'* section.

2.4 Doppler Communications

Select '*Setup*' menu, and then '*Doppler*'.

The following window will appear.

Ensure an Interface Cable is connected. Connect a probe to the Doppler and switch on.

After 2-3 seconds, a message will appear, [LINKUP OK]. The probe frequency, software version and battery voltage will be displayed. If they are not shown, select another Communications Port using the drop down window, and click on '*Linkup*'.

Repeat with other ports until linkup is successful. Click on *OK*.

If the USB-serial adaptor (ACC190) is being used, then select the appropriate 'COM PORT' number.

If the COM PORT number is not known, then:

- Select '*Control Panel*' from the '*START*' menu
- Select '*SYSTEM*'
- Select '*HARDWARE*' Tab
- Select '*DEVICE MANAGER*'
- Select '*PORTS*' (COM & LPT)

The COM PORT number should be shown next to the 'USB-to-serial Com Port (comxx)'.

2.5 Doppler Application

The main window has up to 12 waveform boxes displayed. The highlighted box can then be enlarged by pressing the **<Enter>** key or point to the box using the mouse, and double click.

A Doppler signal can be obtained and displayed by pressing the space bar. To keep the waveform, press the space bar again, or select **'KEEP'** with the mouse.


Press **<Enter>** or click on **OK** to return to the main window. The waveform is shrunk into one of the 12 waveform boxes. The next box in the sequence is automatically highlighted and can be enlarged as before, by pressing **<Enter>**, or point to the box using the mouse and double click.


The two buttons, **MODE** and **START/STOP**, operate in the same way as the **OK** and **KEEP** buttons respectively. This allows the user to obtain signals from various vessels and save them in sequence without having to operate the computer keyboard. The operation of the **MD2/RD2/MD200** with the software, are explained in more detail in the **'Doppler Data Input - Doppler Data Entry'** box in the **'HELP'** facility.

MD2/RD2


MD200


2.6 Doppler Calculations

Various parameters can be calculated on the Doppler waveforms. These are enabled by checking the relevant boxes in the *'Calculations'* window from the *'Setup'* menu.

All parameters are averaged over the number of cardiac cycles shown. For correct calculations, the vertical waveform markers must align with the start of each cardiac cycle. If they do not, then recapture a new display of arterial waveforms. **DR4** will also calculate these parameters on previously stored data obtained from **DR2/DR3**.

If the waveforms below the zero line are not to be included in the calculations, select +ve calc .

2.6.1 Auto Gain Control

If the  is selected, the Doppler waveforms will automatically be adjusted in height after every screen. The waveform height is therefore optimised to between half and full scale.

2.7 VPPG Application

When connected to the RD2, DR4 can display and printout VPPG waveforms.

Select 'VPPG' from the 'Test' menu and attach the VPPG sensor to the leg by following the instructions in the RD2 user manual.


Switch the RD2 on and select the relevant tourniquet box by double clicking or use <Enter>. To start the test, select 'RUN' When the trace has stabilised, (select Zero Draw to observe), a 5 second countdown will begin.

The patient's foot should move in time with the graphic to complete 10 dorsiflexions.

The patient should then rest for 45 seconds and the refill curve is displayed.


The parameter **RT**, **VP** and **HAT** are automatically calculated, but can be altered by selecting the cursors.

The test can be repeated by selecting '**RUN**' again. The traces will change colour, the black trace being active.

The black trace will always be saved after selecting OK.

To complete the VPPG examination, store the curves in the relevant boxes, using the tourniquet cuffs where appropriate. Refer to the RD2 user manual for correct use of the tourniquet cuffs.


- To display the nomogram, select 'NOMOGRAM' from the patient menu. This categorises the results into NORMAL, BORDERLINE and ABNORMAL. (Enable the nomogram from the 'Setup' menu - 'Print Layout').

2.8 ABI Application


When connected to the dopplex Ability, DR4 can display, archive, and printout ABIs and PVR waveforms.

Connect the USB cable supplied with the dopplex Ability, to the unit and the PC. Switch on the Ability and DR4 will automatically find the correct communication port.

Start an ABI test and a new window will be shown.


When the test has finished, data and waveforms will be automatically transferred to the PC.

Select  to save the results or  to 'Save As' the results as a different file name.

Select  to display test history.

Note: the Ability data and history appear as new windows over the Doppler traces. If this window disappears, click on the lower toolbar to allow it to reappear.

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If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor.
In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

Manufactured in the UK by Huntleigh Healthcare Ltd on behalf of;

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden

Huntleigh Healthcare Ltd.

35 Portmanmoor Road, Cardiff, CF24 5HN, United Kingdom
T: +44 (0)29 20485885 sales@huntleigh-diagnostics.co.uk
www.huntleigh-diagnostics.com

Distributed in the USA by:

ArjoHuntleigh, Inc.
2349 West Lake Street, Suite 250
Addison, IL 60101
T: 800-323-1245
www.huntleigh-healthcare.us

1001057-2

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HNE Huntleigh
Industriering Ost 66
47906 Kempen
Germany
T: 02152 551110
F: 02152 551120

1001058-1

Registered No: 942245 England & Wales. Registered Office:
ArjoHuntleigh House, Houghton Hall Business Park, Houghton Regis, Bedfordshire, LU5 5XF
©Huntleigh Healthcare Limited 2005-2014

A Member of the Arjo Family

® and ™ are trademarks of Huntleigh Technology Limited

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