

HUNTLEIGH

A MEMBER OF THE ARJO FAMILY



Smartsigns

Patient Monitoring Accessories
& Consumables Catalogue

2023



CONTENTS

MINIPULSE MP1 **A**

Accessories **3**

Consumables **4**

VITAL SIGNS & PATIENT MONITORS **B**

SpO2 - Own Brand **7**

SpO2 - Nellcor Oximax **8**

NiBP - Reusable Cuffs **9**

NiBP - Single Use Cuffs **11**

ECG **12**

Temperature **14**

Invasive Blood Pressure **16**

TcCO2 **17**

Hardware & Consumables **20**

SUPERSEDED ACCESSORIES **C**

Temperature **24**

BP/Cuffs **26**

SpO2 **27**

ECG/Respiration **30**

Invasive Blood Pressure **31**

CO2 **33**

Hardware & Consumables **37**



ACC-VSM-161

SpO2 Extension Cable
(5'')



ACC-VSM-162

Ear Sensor



ACC-VSM-169

Universal Y Sensor



ACC-VSM-170

Paediatric Finger Sensor
(10-40kg)



ACC-VSM-171M

Reusable SpO2 finger Sensor
Adult (>40kg)



Non-Compatible



Compatible

A

ACC-VSM-163

Adult Finger Sensor
- Box of 10 (>45kg)

**ACC-VSM-164**

Paediatric Finger Sensor
- Box of 10 (15-45kg)

**ACC-VSM-165**

Neonatal Finger Sensor
- Box of 10 (<3kg)

**ACC-VSM-166**

Infant Finger Sensor
- Box of 10 (3-15kg)

**ACC-VSM-167**

Infant Wrap Sensor

**ACC-VSM-168**

Neonatal Wrap Sensor



ACC47

IV Pole Attachment



ACC-VSM-172

Carry Case



ACC-VSM-173

MP1 Desk Stand



ACC-VSM-174

MP1R Desk Stand with built in Charger



SP-750071

Power Adaptor - Specify Country:

MP1R-ADAPT-UK UK

MP1R-ADAPT-EUR Europe

MP1R-ADAPT-US/J USA / Japan

MP1R-ADAPT-ROW Rest of World



ACC-VSM-232

Neonate Tape Attachment
(Pack of 50)

ACC-VSM-233

Infant Tape Attachment
(Pack of 100)



For use with ACC-VSM-167/168



Non-Compatible



Compatible

SpO2 SENSORS

ACC-VSM-289

Adult Re-usable SpO2 Sensor
(10ft)



ACC-VSM-290

SpO2 extension cable
(8ft) - compatible with
ACC-VSM-291
ACC-VSM-292



ACC-VSM-291

Adult Reusable SpO2 sensor - DB9F
(3ft) - compatible with ACC-VSM-290



ACC-VSM-292

Neonate / Paediatric Wrap SpO2 sensor - DB9F
(3ft) - compatible with ACC-VSM-290



ACC-VSM-305

Adult Soft Tip Sensor



B



Non-Compatible



Compatible

ACC-VSM-306

Paediatric Soft Tip Sensor

**ACC-VSM-307**

Infant Soft Tip Sensor

**ACC-VSM-308**

Neonatal Silicone Wrap Sensor

**ACC-VSM-309**DOC10 Compatible Extension cable
(Oximax compatible)**ACC-VSM-310**DS100A Compatible Sensor
(Oximax Compatible)**ACC-VSM-311**Neonate Disposable SpO2 Sensor
(Oximax compatible)

ACC-VSM-312

Paediatric Disposable SpO2 Sensor
(Oximax compatible)

SC300
SC500
SC1200
SC1500

**ACC-VSM-313**

Adult Disposable SpO2 Sensor
(Oximax compatible)

**ACC-VSM-314**

Adult Soft Tip SpO2 Sensor
(Oximax compatible)

**ACC-VSM-315**

Paediatric Soft Tip SpO2 Sensor
(Oximax compatible)

**ACC-VSM-316**

Infant Soft Tip SpO2 Sensor
(Oximax compatible)

**ACC-VSM-317**

Neonatal Silicone Wrap Sensor
(Oximax compatible)






NiBP Cuffs

ACC-VSM-273

Infant Reusable NiBP Cuff
(6-11cm)



ACC-VSM-274

Infant Reusable NiBP Cuff
(10-19cm)



ACC-VSM-275




Paediatric Reusable NiBP Cuff
(18-26cm)



ACC-VSM-276

Adult Reusable NiBP Cuff
(20-28cm)



ACC-VSM-277Adult Reusable NiBP Cuff
(25-35cm)**ACC-VSM-278**Adult Reusable NiBP Cuff
(33-47cm)**ACC-VSM-279**Adult Thigh Reusable NiBP Cuff
(46-66cm)**ACC-VSM-284**NiBP Hose
(3M)**ACC-VSM-280**

Disposable #1 Neonatal Cuff (3-5cm)

**ACC-VSM-281**

Disposable #2 Neonatal Cuff (4-7.6cm)

**ACC-VSM-282**

Disposable #3 Neonatal Cuff (5.6-10.6cm)

ACC-VSM-283

Disposable #4 Neonatal Cuff (7-12.8cm)



ACC-VSM-318

Adult Single Use Cuff
- 5 pcs (45-56cm)

ACC-VSM-319

Adult Single Use Cuff
- 5 pcs (36.5-46cm)

ACC-VSM-320

Adult Single Use Cuff
- 5 pcs (27.5 - 36.5cm)

ACC-VSM-321

Adult Single Use Cuff
- 5 pcs (20.5 - 28.5cm)

ACC-VSM-322

Paediatric Single Use Cuff
- 5 pcs (13.8 - 21.5cm)

ACC-VSM-323

Infant Single Use Cuff
- 5 pcs (9 - 14.8cm)



Non-Compatible




Compatible



ACC-VSM-201

Pre-Wired Neonatal Electrodes
Box of 50 (used with ACC-VSM-303)

SC300	SC500	SC1200	SC1500
☐	☒	☒	☒



ACC-OBS-027

ECG Skin Electrodes
(Box of 300)

☐	☒	☒	☒
--------------------------	-------------------------------------	-------------------------------------	-------------------------------------



ACC-VSM-300

3 Way ECG Cable (IEC)

☐	☒	☒	☒
--------------------------	-------------------------------------	-------------------------------------	-------------------------------------



ACC-VSM-301

5-Way ECG cable (IEC)

SC300	SC500	SC1200	SC1500
☐	☐	☒	☒



ACC-VSM-302

12-Way ECG cable (IEC)

SC300	SC500	SC1200	SC1500
☐	☐	☒	☒



ACC-VSM-303

3-Way Neonatal Trunk

SC300	SC500	SC1200	SC1500
☐	☒	☒	☒



TEMPERATURE



SP-781405

Temperature Probe

SC300	SC500	SC1200	SC1500
☒	☒	☐	☐



SP-781398

Temperature Probe holder

SC300	SC500	SC1200	SC1500
☒	☒	☐	☐



ACC-VSM-324

Adult Surface Temperature Sensor (10K)

SC300	SC500	SC1200	SC1500
☐	☐	☒	☒



ACC-VSM-325

Adult Internal Temperature Sensor (10K)

SC300	SC500	SC1200	SC1500
☐	☐	☒	☒





	SC300	SC500	SC1200	SC1500
ACC-VSM-326 Neonatal Internal Temperature Sensor (10K)	☐	☐	☒	☒
ACC-VSM-327 Neonatal Surface Temperature Sensor (10K)	☐	☒	☒	☒

ACC-VSM-326
Neonatal Internal Temperature Sensor
(10K)

ACC-VSM-327
Neonatal Surface Temperature Sensor
(10K)

Invasive Blood Pressure

ACC-VSM-328

IBP Interface Cable with Abbott connector

ACC-VSM-330

IBP Interface Cable with B&D connector

ACC-VSM-332

IBP Interface Cable with Braun connector

ACC-VSM-334

IBP Interface Cable with Edwards connector

ACC-VSM-336

IBP Interface Cable with Utah connector

ACC-VSM-350

IBP Interface Cable with Medex Logic connector

ACC-VSM-329

Disposable Pressure transducer Kit with Abbott connector

ACC-VSM-331

Disposable Pressure transducer Kit with B&D connector

ACC-VSM-333

Disposable Pressure transducer Kit with Braun connector

ACC-VSM-335

Disposable Pressure transducer Kit with Edwards connector

ACC-VSM-337

Disposable Pressure transducer Kit with Utah connector

ACC-VSM-350

Disposable Pressure transducer Kit with Medex Logic connector

SC300
SC500
SC1200
SC1500



EtCO₂**ACC-VSM-203**

Adult Nasal Cannula (pack of 10)

ACC-VSM-204

Paediatric Nasal Cannula (pack of 10)

ACC-VSM-205

Infant / Neonatal Nasal Cannula (pack of 10)

ACC-VSM-338Adult Nasal Cannula with filter - High humidity
(Pack of 10)**ACC-VSM-339**Adult Nasal Cannula with filter - High humidity
(Pack of 10)**ACC-VSM-340**Adult Nasal Cannula with filter - High humidity
(Pack of 10)**ACC-VSM-206**Straight Sampling line
- Side stream (pack of 10)SC300
SC500
SC1200
SC1500

Non-Compatible



Compatible

ACC-VSM-207

Straight Sampling line with filter
- Side Stream (Pack of 10)

**ACC-VSM-208**

Adult / Paediatric Airway Adaptor Kit -
Side Stream (pack of 10)

**ACC-VSM-210**

Adult / Paediatric Airway Adaptor Kit
with filter
- Side Stream (pack of 10)

**ACC-VSM-209**

Infant / Neonatal Airway Adaptor Kit
- Side Stream (Pack of 10)

**ACC-VSM-211**

Infant / Neonatal Airway Adaptor Kit
with filter
- Side Stream (pack of 10)

**ACC-VSM-341**

CAPNOSTAT 5 Main Stream Sensor

**ACC-VSM-342**

Main stream Sensor Interface cable

**ACC-VSM-343**

Disposable Airway Adapter
- Main Stream Adult/Paediatric (pack of 10)



ACC-VSM-344

Disposable Airway Adapter
- Main Stream Infant / Paediatric (pack of 10)

ACC-VSM-345

Reusable Airway Adapter
- Main Stream Adult/Paediatric

ACC-VSM-346

Reusable Airway Adapter
- Main Stream Infant / Paediatric

SC300
SC500
SC1200
SC1500



Hardware & Consumables



ACC-VSM-153

Light Duty Roll Stand with Integrated Storage Basket - Requires Fixing kit

SC300
SC500
SC1200
SC1500



ACC-VSM-154

Wall Mount - Requires Fixing kit



ACC-VSM-187

Utility Hook - Requires ACC VSM 154



ACC-VSM-189

3.5" Utility Basket - Requires ACC VSM 154



ACC-VSM-294


IV Pole Clamp



ACC-VSM-288

Fixing Kit

- Use with ACC VSM 153 & 154

**ACC-VSM-348**

Earth Cable

**ACC-VSM-353**

Fixing Kit

- Used with ACC VSM 153 & 154

**ACC-VSM-349**

Fixing Kit

- Used with ACC VSM 153 & 154

**ACC-VSM-286**

Lens Filter for Thermometer

- Box of 40 (800pcs)

ACC-VSM-287

Lens Filter for Thermometer

- Box of 400 (8,000pcs)

ACC-VSM-293

Lens Filter for Thermometer

- Box of 10 (200pcs)

**ACC-VSM-347**

Printer Paper with grid

- 50mm (20m)




LEGACY ACCESSORIES



Temperature

ACC-VSM-43

Oral Probe - single (MK1 module only)



ACC-VSM-44

Rectal Probe - single (MK1 module only)



ACC-VSM-193

Suretemp Plus Oral probe - Single (V3)



ACC-VSM-194

Suretemp Plus Rectal probe - Single (V3)



ACC-VSM-11

Central - Adult

ACC-VSM-12

Central - Infant

ACC-VSM-13

Central - Neonate



C

☐ Non-Compatible

☒ Compatible

ACC-VSM-14

Surface - Adult

ACC-VSM-15

Surface - Infant

ACC-VSM-16

Surface - Neonate

BP / Cuffs

ACC-VSM-223

NiBP Cuff - Infant

Box of 5 (8-13cm)

ACC-VSM-224

NiBP Cuff - Child

Box of 5 (12-19cm)

ACC-VSM-225

NiBP Cuff - Small Adult

Box of 5 (17-25cm)

ACC-VSM-129

NiBP Cuff - Neonatal

Box of 10 (Size 1: 3-6cm)

ACC-VSM-130

NiBP Cuff - Neonatal

Box of 10 (Size 2: 4-8cm)

ACC-VSM-132

NiBP Cuff - Neonatal

Box of 10 (Size 4: 7-13cm)

ACC-VSM-133

NiBP Cuff - Neonatal

Box of 10 (Size 5: 8-15cm)

ACC-VSM-195

All Purpose Cuffs - Infant (8-13cm)

ACC-VSM-196

All Purpose Cuffs - Child (12-19cm)

ACC-VSM-197

All Purpose Cuffs - Small Adult (17-25cm)

ACC-VSM-198

All Purpose Cuffs - Adult (23-33cm)

ACC-VSM-199

All Purpose Cuffs - Large Adult (31-40cm)

ACC-VSM-200

All Purpose Cuffs - Thigh Adult (38-50cm)

SL600
SC750
SC1000



Non-Compatible



Compatible

BP Hoses

ACC-VSM-142

NiBP Hose - 6ft Adult / Paediatric
(MK 1 product)

ACC-VSM-143

NiBP Hose - 8ft Adult / Paediatric
(MK 1 product)

ACC-VSM-144

NiBP Hose - 8ft Neonatal
(MK 1 product)

ACC-VSM-184

Neonatal Hose 3.5m
(MK 2 & MK 3 product)

ACC-VSM-185

Adult Hose 3.5m
(MK 2 & MK 3 product)

ACC-VSM-229

Hose for Disposable Cuffs (ACC VSM 223-228)
(MK 2 & MK 3 product)

SL600
SC750
SC1000



ACC-NEL-32203231

Nellcor OxiMAX dura-Y Multi-site Sensor
Adult (1m / 3ft Cable)

SL600	SC750	SC1000
☒	☐	☐



ACC-MAX-A (D-25)

Disposable OxiMAX SpO2 sensor
Adult (box of 24)

☒	☐	☐
-------------------------------------	--------------------------	--------------------------



ACC-MAX-I (I-20)

Disposable OxiMAX SpO2 sensor
Infant (box of 24)

☒	☐	☐
-------------------------------------	--------------------------	--------------------------



ACC-MAX-N (N-25)

Disposable OxiMAX SpO2 sensor
Neonatal (box of 24)

☒	☐	☐
-------------------------------------	--------------------------	--------------------------



ACC-MAX-P (D-20)

Disposable OxiMAX SpO2 sensor
Paediatric (box of 24)

☒	☐	☐
-------------------------------------	--------------------------	--------------------------



ACC-VSM-162

Reusable Ear Sensor

**ACC-VSM-167**

Reusable Infant Wrap Sensor

**ACC-VSM-168**

Reusable Neonatal Wrap Sensor

**ACC-VSM-169**

Reusable Universal Y Sensor

**ACC-VSM-170**Reusable Paediatric Finger Sensor
(5-15kg)**ACC-VSM-171**

Reusable Comfort Clip Finger Sensor



ACC-VSM-188

BCI Extension Cable - 8"

SL600	SC750	SC1000
☐	☒	☒



ACC-VSM-163

BCI Disposable Adult Finger Sensor
Box of 10 (>45kg)

☐	☒	☒
--------------------------	-------------------------------------	-------------------------------------



ACC-VSM-164

BCI Disposable Paediatric Finger Sensor
Box of 10 (15-45kg)

☐	☒	☒
--------------------------	-------------------------------------	-------------------------------------



ACC-VSM-165

BCI Disposable Neonatal Finger Sensor
Box of 10 (<3kg)

☐	☒	☒
--------------------------	-------------------------------------	-------------------------------------



ACC-VSM-166

BCI Disposable Infant Finger Sensor
Box of 10 (3-15kg)

☐	☒	☒
--------------------------	-------------------------------------	-------------------------------------



ECG / Respiration

ACC-VSM-01
3- Way Trunk

SL600	SC750	SC1000
☐	☒	☒



ACC-VSM-05
3 Way Combined ECG / RESP patient cable

☐	☒	☒
--------------------------	-------------------------------------	-------------------------------------



ACC-VSM-06
3 Way Fly Lead set RED/GREEN/YELLOW

☐	☒	☒
--------------------------	-------------------------------------	-------------------------------------



ACC-VSM-146
5 Way Trunk

☐	☒	☒
--------------------------	-------------------------------------	-------------------------------------



ACC-VSM-147
5 Way Combined ECG/RESP patient cable

☐	☒	☒
--------------------------	-------------------------------------	-------------------------------------



ACC-VSM-148
5 Way Fly lead Set


☐	☒	☒
--------------------------	-------------------------------------	-------------------------------------



ACC-VSM-201

Pre wired Neonatal Electrodes
(Box of 50)

SL600	SC750	SC1000
☐	☒	☒



Invasive Blood Pressure

ACC-VSM-26

Re-usable Transducer

SL600	SC750	SC1000
☐	☐	☒



ACC-VSM-27

Single Line Monitoring Kit
RED (pack of 20)

☐	☐	☒
--------------------------	--------------------------	-------------------------------------



ACC-VSM-28

Single Line Monitoring Kit
BLUE (pack of 20)

☐	☐	☒
--------------------------	--------------------------	-------------------------------------



ACC-VSM-029

Dual Line Monitoring Kit
(Pack of 20)

SL600
SC750
SC1000

**ACC-VSM-30**

Blood Pressure Cable

**ACC-VSM-31**

Mounting Bracket for Two Transducers
(Used with ACC-VSM-32)

**ACC-VSM-32**

Stand Clamp
(Used with ACC-VSM-31)

**ACC-VSM-86**

Single Line IBP Cable
- Edwards



ACC-VSM-87

Single Line IBP Cable
- B&D

SL600
SC750
SC1000

**ACC-VSM-104**

Single Line IBP Cable
- Abbott

ACC-VSM-109

Single Line IBP Cable
- MEDEX (MX9504 series)

ACC-VSM-110


Single Line IBP Cable
- KIMAL / ARROW

**ACC-VSM-221**

IBP Monitor Cable
- Tanspac

ACC-VSM-222

IBP Monitor Cable
- COMBI Trans

**CO2 - LOFLO SIDESTREAM****ACC-VSM-203**

Respironics Adult Nasal Cannula
(Box of 10)

ACC-VSM-204

Respironics Paediatric Nasal Cannula
(Box of 10)

ACC-VSM-205

Respironics Infant / Neonatal Nasal Cannula
(Box of 10)

SL600
SC750
SC1000



Non-Compatible



Compatible

ACC-VSM-206

Respironics Straight Sampling Line
(Pack of 10)

ACC-VSM-207

Respironics Straight Sampling Line with Filter
(Pack of 10)

ACC-VSM-208

Respironics Airway Adaptor Kit
Adult / Paediatric >4mm (pack of 10)

ACC-VSM-210

Respironics Airway Adaptor Kit with Filter
Adult / Paediatric >4mm (pack of 10)

ACC-VSM-209

Respironics Airway Adaptor Kit
Infant / Neonatal <4mm (pack of 10)

ACC-VSM-211

Respironics Airway Adaptor Kit with Filter
Infant / Neonatal <4mm (pack of 10)


SL600
SC750
SC1000

**CO2 -Calibration Only - Welsh Allyn Techno**


SL600
SC750
SC1000

**ACC-VSM-33**

CO2 Absorber for Calibration
(Pack of 5)

**ACC-VSM-34**

Water Trap
(Pack of 10)



SURE CO2 - Welsh Allyn Technology

ACC-VSM-36

EtCO2 Nasal Sampline Line (soft)
Adult (box of 25)

SL600
SC750
SC1000



ACC-VSM-38

90o Elbow
(Box of 50)



ACC-VSM-96

EtCO2 Dryer Sampling Line
HH (box of 25)

ACC-VSM-97

EtCO2 Dryer Sampling Line
HH with TEE (box of 25)

ACC-VSM-98

EtCO2 Dryer Sampling Line
LH (box of 25)

ACC-VSM-99

EtCO2 Dryer Sampling Line
LH with TEE (box of 25)



ACC-VSM-100

EtCO2 Nasal Sampling Line (Soft)
Paediatric (box of 25)



ACC-VSM-101

EtCO2 Nasal Sampling Line (Soft)
Infant (box of 25)



Non-Compatible



Compatible

ACC-VSM-115

Adult / Paediatric TEE with Sampling Port
(Pack of 50)

ACC-VSM-116

Infant connector with side port
2.5mm (box of 50)

ACC-VSM-117

Infant connector with side port
3.0mm (box of 50)

ACC-VSM-118

Infant connector with side port
3.5mm (box of 50)

ACC-VSM-119

Infant connector with side port
4.0mm (box of 50)



CO2 - Delivery & Sampling Lines

ACC-VSM-35

O2/EtCO2 Nasal Sampling Line (Soft)
Adult (box of 25)



ACC-VSM-102

O2/EtCO2 Nasal Sampling Line (Soft)
Paediatric (box of 25)



ACC-VSM-103

O2/EtCO2 Nasal Sampline Line (Soft)
Infant (box of 25)



Hardware

ACC-VSM-04

Soft Carry Bag



ACC-VSM-95

Carry Bag



ACC-VSM-151

Fixing kit

(Used with ACC VSM153 & 154)



ACC-VSM-152

Fixing kit

(Used with ACC VSM 153 & 154)



ACC-VSM-153

Light Duty Roll Stand

(With integrated storage basket)



ACC-VSM-154

Wall Mount (complete kit)



ACC-VSM-159

Fixing kit
(Used with ACC VSM 153 & 154)



ACC-VSM-160

DIN Rail Mount



ACC-VSM-175

Horizontal Rail Mount



ACC-VSM-187

Utility Hook
(Compatible with ACC VSM 154)



ACC-VSM-189

3.5" Utility Basket
(Compatible with ACC VSM 154)



Consumable

ACC179
Adhesive Strap
(Pack of 240)

ACC-VSM-45
Probes Covers
(Pack of 1000)


ACC-VSM-39
Chart Paper
(Pack of 10)

ACC-VSM-40
Chart Paper
(Pack of 100)

ACC-VSM-145
Thermal Paper
SL600 units(pack of 3)

ACC-VSM-232
Neonatal Tape Attachment
(Pack of 50)

ACC-VSM-233
Infant Tape Attachment
(Pack of 50)



As a proud member of the Arjo family, we have been committed to supporting healthcare professionals in improving outcomes and enhancing patient wellbeing since 1979. We do this through our proven solutions for Vascular Assessment & Treatment and Fetal & Patient Monitoring. With innovation and customer satisfaction as our guiding principles, we strive for clinical excellence and improved performance, for life.

Manufactured and distributed by Huntleigh Healthcare Ltd.

35 Portmanmoor Road, Cardiff, CF24 5HN, United Kingdom
T: +44 (0)29 20485885 sales@huntleigh-diagnostics.co.uk
www.huntleigh-diagnostics.com

Registered No: 942245 England & Wales. Registered Office:
ArjoHuntleigh House, Houghton Hall Business Park, Houghton Regis, Bedfordshire, LU5 5XF
©Huntleigh Healthcare Limited 2019

A Member of the Arjo Family

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

HUNTLEIGH

