

HUNTLEIGH

A MEMBER OF THE ARJO FAMILY



Sonicaid

Fetal Monitoring Accessories &
Consumables Catalogue

2024



Contents

<u>HANDHELD DOPPLERS</u>	A
Accessories	2
Consumables	3
<u>TEAM3</u>	B
Accessories	4
Hardware	7
Software	8
Compatibility	15
<u>WIRELESS</u>	C
Accessories	12
Consumables	14
Compatibility	15
<u>LEGACY ACCESSORIES</u>	D
Dopplers	17
Monitors Accessories	19
Monitors Hardware	24
Monitors Consumables	26





DOPPLERS

ACC21
Stereo Headset

SR2	SR3	SRX	D920	D930	FD1	FD3
Compatible	Compatible	Compatible	Compatible	Compatible	Compatible	Compatible



ACC34
Soft Carry Pouch

Compatible	Compatible	Compatible	Compatible	Compatible	Compatible	Compatible
------------	------------	------------	------------	------------	------------	------------



ACC226
Medical Grade Rechargeable Kit

Compatible	Compatible	Compatible	Non-Compatible	Non-Compatible	Non-Compatible	Non-Compatible
------------	------------	------------	----------------	----------------	----------------	----------------



ACC232
Dopplex Stand - shipped flat packed Includes adaptor plate

Compatible	Compatible	Compatible	Compatible	Compatible	Compatible	Compatible
------------	------------	------------	------------	------------	------------	------------



Compatible
 Non-Compatible

ACC227

Micro SD Card

SR2	SR3	SRX	D920	D930	FD1	FD3
Compatible	Compatible	Compatible	Non-Compatible	Non-Compatible	Non-Compatible	Non-Compatible



ACC-OBS-080

Splash-proof (IPX2) Doppler Pouch

Compatible	Compatible	Compatible	Non-Compatible	Non-Compatible	Non-Compatible	Non-Compatible
------------	------------	------------	----------------	----------------	----------------	----------------



ACC3

Ultrasonic Gel
(Box of 12 x 250ml Bottles)

Compatible	Compatible	Compatible	Compatible	Compatible	Compatible	Compatible
------------	------------	------------	------------	------------	------------	------------



ACC24

Ultrasonic Gel
(Box of 12 x 60ml Tubes)

Compatible	Compatible	Compatible	Compatible	Compatible	Compatible	Compatible
------------	------------	------------	------------	------------	------------	------------



ACC-1300-0154

Ultrasonic Gel
(5lt container)

Compatible	Compatible	Compatible	Compatible	Compatible	Compatible	Compatible
------------	------------	------------	------------	------------	------------	------------





TEAM3

DOPPLERS

TEAM3

WIRELESS

LEGACY

ACC-OBS-008

1.0MHz Ultrasound Transducer
Green (2 meters)

ACC-OBS-008L

1.0MHz Ultrasound Transducer
Green (3 meters)

ACC-OBS-009

Toco Transducer
Pink (2 meters)

ACC-OBS-009L

Toco Transducer
Pink (3 meters)

SP 7775-6901

Event Marker

ACC-OBS-066

FECG Leg Plate
Safelinc

ACC-OBS-068

FECG Leg Plate

Intrapartum
Antepartum



Compatible

Non -Compatible

B

ACC-OBS-069

FECG Leg Plate
Philips Compatible

Intrapartum
Antepartum



ACC-OBS-120

Extra Large Conical BP Cuff
(38-46cm Circumference)



ACC-OBS-121

Large BP Cuff
(32-42cm Circumference)



ACC-OBS-122

Medium BP Cuff
(24-32cm Circumference)



ACC-OBS-123

Small BP Cuff
(17-26cm Circumference)



ACC-OBS-124

Hose for BP Cuffs



ACC-OBS-071

SpO2 Cable
(Compatible with BCI sensors)



ACC-OBS-126

Comfort Clip Flinger Sensor
BCI compatible

Intrapartum
Antepartum



ACC-NEL-DS-100A

SpO2 Sensor
Nellcor OxiMax
(Compatible with ACC-OBS-071)



ACC-NEL-DOC-10

SpO2 Interface Cable
Nellcor OxiMax (3m / 10Ft Cable)



ACC-OBS-070

Maternal ECG Lead



ACC-8400- 6937

Intra Plus Interconnection Lead



ACC-OBS-032

Koala Interconnection Lead



HARDWARE

ACC-OBS-072

Trolley / Cart with fixing kit
(Requires ACC-OBS-083)

Intrapartum
Antepartum




ACC-OBS-076

Universal Wall Mount
(Requires ACC-OBS-090)



ACC-OBS-090

Universal Wall Mount Fixing plate
(Used with ACC-OBS-076)



ACC-OBS-077

Paper Tray - Philips

ACC-OBS-078

Paper Tray - Corometrics (GE)

ACC-OBS-079

Paper Tray - Huntleigh



ACC-OBS-083

Fixing plate
(Used with ACC-OBS-072)



ACC-OBS-091

Carry bag



CONSUMABLES

ACC3

Ultrasonic Gel
(Box of 12 x 250ml bottles)

Intrapartum
Antepartum



ACC24

Ultrasonic Gel
(Box of 12 x 60ml Tubes)



ACC-MI1136

Reusable Transducer belts (Latex Free)
Pair



ACC220

Disposable Transducer belts (Latex Free)
Pack of 100 pairs



ACC-321414

ArchiTrace 25 year archival Chart Paper
45m Length (pack of 20)



ACC-OBS-109

FECG Scalp Electrodes
Neoventa connect (box of 50)



Compatible


Non -Compatible

ACC-OBS-027

Maternal ECG Skin Electrodes
(5mm Stud - box of 300)

Intrapartum
Antepartum


☒	☐
-------------------------------------	--------------------------



ACC-7000AA0

Spiral Electrodes
QWIK connect (box of 50)


☒	☐
-------------------------------------	--------------------------



ACC-OBS-020

FECG Copeland Electrodes
QWIK connect (box of 25)

☒	☐
-------------------------------------	--------------------------



ACC-2464AA0

Adhesive Leg Electrodes
- QWIK Connect (3mm Stud - box of 50)

☒	☐
-------------------------------------	--------------------------



ACC-900X259

FECG Spiral Electrodes
Safelinc (box of 50)


☒	☐
-------------------------------------	--------------------------



ACC-900X260

Reference Electrodes
Safelinc (box of 50)

☒	☐
-------------------------------------	--------------------------



ACC-8400-8011
Intran Plus IUP Catheter
(Box of 10)

ACC-1300-145
Ultrasonic Gel
20g Sachet (box of 48)

ACC-1300-0154
Ultrasonic Gel
5lt Container (each)

Intrapartum
Antepartum



Compatible



Non -Compatible

WIRELESS



ACC-OBS-093

Interface Cable
(Wireless to TEAM3)



ACC-OBS-095

Ultrasound Transducer
(European RF Standard)



ACC-OBS-096

Toco-T Transducer
(European RF Standard)



ACC-OBS-099

Toco-E Transducer with MEGG/
FECG Connector
(European RF Standard)



ACC-OBS-101

2-lead MEGG cable
(IEC) (requires ACC-OBS-027)



C




Compatible



Non -Compatible

ACC-OBS-102


FECG Connection Cable
(Tyco Safelinc)

**ACC-OBS-106**

FECG Connection Cable
(Philips/Neovanta)

**ACC-OBS-107**

FECG Connection Cable
(GE QWIK)

**ACC-OBS-027**

Maternal ECG Skin Electrodes
- 5mm Stud (box of 300) (used with ACC-OBS-101)

**ACC-900X259**

FECG Spiral Electrodes
Safelinc (box of 50)

**ACC-700AAA0**


FECG Scalp Electrodes
Qwik connect (box of 50)



ACC-OBS-020


FECG Scalp Electrodes

Rocket Copeland connect (box of 50)

**ACC-OBS-109**

FECG Scalp Electrodes

Neoventa connect (box of 50)

**ACC3**

Ultrasonic Gel

(Box of 12 x 250ml Bottles)

**ACC24**

Ultrasonic Gel

(Box of 12 x 60ml Tubes)

**ACC-1300-145**

Ultrasonic Gel

20g Sachet (box of 48)

**ACC-1300-0154**

Ultrasonic Gel

5lt Container (each)



ACC-MI1136

Reusable Transducer belts (Latex Free)

Pair

**ACC220**

Disposable Transducer belts (Latex Free)

Pack of 100 pairs

**ACC-OBS-108**

Mains Y Splitter Power Cable

(1m)

**ACC-OBS-103**

Trolley Bracket kit

(Used with ACC-OBS-072)

**ACC-OBS-114**

FTS-3 Wall Mount

(Used with ACC-OBS-076)







COMPATIBILITY

LEGACY

WIRELESS

TEAM3

DOPPLERS





ACC900X259
SAFELINC FECG ELECTRODE

ACC900X260
SAFELINC REFERENCE ELECTRODE



ACC-700AAO
FSE - QWIK CONNECT SPIRAL

ACC-2464AAO
QWIK CONNECT REFERENCE ELECTRODES



ACC-OBS-068
FECG/DECG QWIK CONNECTION LEAD PLATE



ACC-OBS-020
FSE - ROCKET COPELAND



■ Compatible

□ Non -Compatible

LEGACY ACCESSORIES



ACC52

Doppler Mobile Stand
(Excluding Sonicaid One)

ACC230

Adaptor Plate

ACC171

Rechargeable Battery Kit (UK)

ACC172

Rechargeable Battery Kit (USA)

ACC173

Rechargeable Battery Kit (EU)

ACC180

Pole Stand
(MD200 only)

046V002

Carry pouch
(Sonicaid One only)

Dopplers
B4000XS
TEAM
FM800 Encore
Freedom



D



Compatible



Non -Compatible

US1

Ultrasound Transducer

US1L

Ultrasound Transducer
(3m Length)

CT1

Contraction Transducer

CT1L

Contraction Transducer
(3m Length)

EM1

Patient Event Marker

EM1L

Patient Event Marker
(3m Length)

ACC30

Carry Case


ACC-MI1136

Reusable Transducer Belts
- Latex Free (Pair)

ACC-739087

FECG Connection Lead
- QWIK connect

Dopplers
B4000XS
TEAM
FM800 Encore
Freedom



8400-6919

Ultrasound Transducer
- 1.5MHz (Yellow)

8400-6919L

Ultrasound Transducer
- 1.5MHz (Yellow) (2.8m Length)

8400-6920

Ultrasound Transducer
- 2MHz (Blue)

8400-6920L

Ultrasound Transducer
- 2MHz (Blue) (2.8m Length)

ACC-OBS-009

Toco Transducer
(Pink)

ACC-OBS-009L

Toco Transducer
(Pink) (3m Length)

ACC-8402-6969

Maternal ECG Lead
(TeamIP Only)

ACC-017D065


FECG Connection Lead
- Safelinc Connect

ACC-NEL-DS-100A

Comfort Clip Finger Sensor
- Nellcor OxiMax Dual Sensor

Dopplers
B4000XS
TEAM
FM800 Encore
Freedom

☐	☐	☒	☐	☐
☐	☐	☒	☐	☐
☐	☐	☒	☐	☐
☐	☐	☒	☐	☐
☐	☐	☒	☐	☐
☐	☐	☒	☐	☐
☐	☐	☒	☐	☐



SP 7775-6901

Event Marker

Dopplers	B4000XS	TEAM	FM800 Encore	Freedom
☐	☐	☒	☒	☐



ACC-8400-6937

Interconnection Lead

- Intra Plus

☐	☐	☒	☒	☐
--------------------------	--------------------------	-------------------------------------	-------------------------------------	--------------------------



ACC-VSM-171

Comfort Clip Finger Sensor

- BCI

☐	☐	☒	☒	☐
--------------------------	--------------------------	-------------------------------------	-------------------------------------	--------------------------



ACC-OBS-008

1.0MHz Ultrasound Transducer
Green (2 meters)

☐	☐	☐	☒	☐
--------------------------	--------------------------	--------------------------	-------------------------------------	--------------------------



ACC-OBS-008L

1.0MHz Ultrasound Transducer
Green (3 meters)

☐	☐	☐	☒	☐
--------------------------	--------------------------	--------------------------	-------------------------------------	--------------------------



ACC-OBS-009

Toco Transducer
Pink (2 meters)

☐	☐	☐	☒	☐
--------------------------	--------------------------	--------------------------	-------------------------------------	--------------------------

ACC-OBS-009L

Toco Transducer
Pink (3 meters)

☐	☐	☐	☒	☐
--------------------------	--------------------------	--------------------------	-------------------------------------	--------------------------



ACC-OBS-007

FECG Leg Plate Adaptor
- Philips Connector

☒ Compatible

☐ Non -Compatible

ACC-OBS-006

SpO2 Cable Connector
(Plastic Connector)

<i>Dopplers</i>	<i>B4000XS</i>	<i>TEAM</i>	<i>FM800 Encore</i>	<i>Freedom</i>
☐	☐	☐	☒	☐



ACC-OBS-021

FECG Connection Lead
- Eziplug

☐	☐	☐	☒	☐
--------------------------	--------------------------	--------------------------	-------------------------------------	--------------------------



ACC-OBS-022

FECG Leg Plate Adaptor
- Safelinc Connector

☐	☐	☐	☒	☐
--------------------------	--------------------------	--------------------------	-------------------------------------	--------------------------



ACC-OBS-023

FECG Leg Plate Adaptor
- QWIK Connector

☐	☐	☐	☒	☐
--------------------------	--------------------------	--------------------------	-------------------------------------	--------------------------



ACC-OBS-024

Maternal ECG Lead

☐	☐	☐	☒	☐
--------------------------	--------------------------	--------------------------	-------------------------------------	--------------------------



ACC-OBS-032

Interconnection Lead
- Koala

☐	☐	☐	☒	☐
--------------------------	--------------------------	--------------------------	-------------------------------------	--------------------------



ACC-OBS-042

SpO2 Cable Connector
(Metal Connector)

<i>Dopplers</i>	<i>B4000XS</i>	<i>TEAM</i>	<i>FM800 Encore</i>	<i>Freedom</i>
☐	☐	☐	☒	☐



ACC-MI1184

Adult Cuff
(31-40cm)
(Requires ACC-MI1160)

☐	☐	☐	☒	☐
--------------------------	--------------------------	--------------------------	-------------------------------------	--------------------------



ACC-MI1160

Hose for Cuffs
(Used with ACC-MI1159 / MI1184)

☐	☐	☐	☒	☐
--------------------------	--------------------------	--------------------------	-------------------------------------	--------------------------



ACC-MI1159

Adult Cuff
(23-33cm)
(Requires ACC-MI1160)

☐	☐	☐	☒	☐
--------------------------	--------------------------	--------------------------	-------------------------------------	--------------------------



ACC-OBS-056

Interface Cable
Connection between FM and
Freedom


☐	☐	☐	☒	☒
--------------------------	--------------------------	--------------------------	-------------------------------------	-------------------------------------



ACC-OBS-050

Mounting System
(Use with FM800 Encore only)


<i>Dopplers</i>	<i>B4000XS</i>	<i>TEAM</i>	<i>FM800 Encore</i>	<i>Freedom</i>
☐	☐	☐	☒	☒



ACC-OBS-051

Transducer Belt Clip
(Pack of 5 pairs)

☐	☐	☐	☐	☒
--------------------------	--------------------------	--------------------------	--------------------------	-------------------------------------



ACC-OBS-052

Ultrasound Transducer

☐	☐	☐	☐	☒
--------------------------	--------------------------	--------------------------	--------------------------	-------------------------------------



ACC-OBS-053

Toco Transducer

☐	☐	☐	☐	☒
--------------------------	--------------------------	--------------------------	--------------------------	-------------------------------------



HARDWARE

ACC-OBS-002

TEAM Trolley with Fixing Kit

ACC-OBS-003

FM Trolley with Fixing Kit

ACC-OBS-004

BD Trolley with Fixing Kit

ACC-OBS-072

Team with Freedom Bracket Trolley
With Fixing Kit

ACC31

Wall Bracket

ACC-OBS-037

Paper Tray - Philips

ACC-OBS-038

Paper Tray - Corometrics (GE)

ACC-OBS-039

Paper Tray - Huntleigh

ACC-OBS-046

UPS Power Supply
(230V)

ACC-OBS-047

UPS Power Supply
(120V)

ACC-OBS-048

Trolley Tray

Dopplers
B4000XS
TEAM
FM800 Encore
Freedom



 **Compatible**

 **Non -Compatible**

Dopplers
B4000XS
TEAM
FM800 Encore
Freedom

ACC-OBS-082
Fixing kit

☐	☐	☐	☐	☒
--------------------------	--------------------------	--------------------------	--------------------------	-------------------------------------



ACC-OBS-084
Aerial Mounting Kit
(Cable, Holder Clip & Allen Key)

☐	☐	☐	☐	☒
--------------------------	--------------------------	--------------------------	--------------------------	-------------------------------------



ACC-OBS-085
Trolley Fixing Kit

☐	☐	☐	☒	☐
--------------------------	--------------------------	--------------------------	-------------------------------------	--------------------------



ACC-OBS-086
Fixing plate
(Used with ACC-OBS-072)

☐	☒	☐	☐	☐
--------------------------	-------------------------------------	--------------------------	--------------------------	--------------------------



ACC223
Universal Wall Mount

☐	☐	☐	☒	☐
--------------------------	--------------------------	--------------------------	-------------------------------------	--------------------------



CONSUMABLES

ACC3

Ultrasonic Gel
(Box 12 x 250ml bottles)

<i>Dopplers</i>	<i>B4000XS</i>	<i>TEAM</i>	<i>FM800 Encore</i>	<i>Freedom</i>
☒	☒	☒	☒	☐



ACC14

Dopplex Printa II Paper
(Pack of 2, 28m per roll)

☒	☐	☐	☐	☐
-------------------------------------	--------------------------	--------------------------	--------------------------	--------------------------



ACC15

Plain Paper
(Box of 10)

☐	☒	☐	☐	☐
--------------------------	-------------------------------------	--------------------------	--------------------------	--------------------------



ACC24

Ultrasonic Gel
(Box of 12 x 60ml Tubes)

☒	☒	☒	☒	☐
-------------------------------------	-------------------------------------	-------------------------------------	-------------------------------------	--------------------------



ACC66

Twins Paper
(Box of 10)

☐	☒	☐	☐	☐
--------------------------	-------------------------------------	--------------------------	--------------------------	--------------------------



ACC220

Disposable Latex Free Transducer Belts
(Pack of 100 wwpairs)

☐	☒	☒	☒	☐
--------------------------	-------------------------------------	-------------------------------------	-------------------------------------	--------------------------




☒ Compatible

☐ Non -Compatible

Dopplers
B4000XS
TEAM
FM800 Encore
Freedom


ACC-OBS-109
FECG Single Helix Spiral Electrodes
- Neovental (box of 75)

☐	☒	☐	☐	☐
--------------------------	-------------------------------------	--------------------------	--------------------------	--------------------------




ACC-OBS-014
FECG Electrodes
- Copeland bare wires (Box of 25)

☐	☐	☐	☒	☐
--------------------------	--------------------------	--------------------------	-------------------------------------	--------------------------



ACC-OBS-020
Copeland Fetal Scalp Electrode
-Qwik connect (box of 25)

☐	☐	☐	☒	☐
--------------------------	--------------------------	--------------------------	-------------------------------------	--------------------------



ACC-OBS-027
Maternal ECG skin electrodes
(5mm Stud - Box of 300)

☐	☐	☒	☒	☐
--------------------------	--------------------------	-------------------------------------	-------------------------------------	--------------------------




ACC-OBS-054
Silicone Sheath
(Each)

☐	☐	☐	☐	☒
--------------------------	--------------------------	--------------------------	--------------------------	-------------------------------------



ACC-1300-0216
Membrane for Toco Transducer
(Roll of 50)

☐	☐	☒	☐	☐
--------------------------	--------------------------	-------------------------------------	--------------------------	--------------------------



ACC-1300-0145

Ultrasonic Gel
(Box of 48 x 20gms Sachet)

Dopplers
B4000XS
TEAM
FM800 Encore
Freedom

**ACC-1300-0154**


Ultrasonic Gel
(5lt container)

**ACC-900X260**

Reference Electrodes
- Safelinc (Box of 50)

**ACC-2464AAO**


Adhesive Leg Electrodes
- Qwik Connect (3mm Stud - box of 500)

**ACC-321414**

ArchiTrace 25 year Archival Chart Paper
(Pack of 20)

**ACC-700AAO**

Spiral electrodes
- Qwik connect (box of 50)



ACC-8400-8011

IUP Catheter
- Intran Plus (Box of 10)

Dopplers
B4000XS
TEAM
FM800 Encore
Freedom



ACC-900X259

FECG Spiral Electrodes
- Safelinc (box of 50)



As a proud member of the Arjo family, we have been committed to supporting healthcare professionals in improving outcomes and enhancing patient wellbeing since 1979. We do this through our proven solutions for Vascular Assessment & Treatment and Fetal & Patient Monitoring. With innovation and customer satisfaction as our guiding principles, we strive for clinical excellence and improved performance, for life.

Manufactured and distributed by Huntleigh Healthcare Ltd.

35 Portmanmoor Road, Cardiff, CF24 5HN, United Kingdom
T: +44 (0)29 20485885 sales@huntleigh-diagnostics.co.uk
www.huntleigh-diagnostics.com

Registered No: 942245 England & Wales. Registered Office:
ArjoHuntleigh House, Houghton Hall Business Park, Houghton Regis, Bedfordshire, LU5 5XF
©Huntleigh Healthcare Limited 2019

A Member of the Arjo Family

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

AW:: 1001043-4